

Product Datasheet

HBcAg(1H8) antibody (Biotin) (orb457769)

Description	HBcAg(1H8) antibody (Biotin)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombined human hepatitis B core antigen protein
Target	HBcAg
Preservatives	0.01M PBS(pH7.4) with 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q89714
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.