

Product Datasheet

E2 tag antibody (Biotin) (orb457888)

| | |
|----------------------------|--|
| Description | E2 tag antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | All |
| Conjugation | Biotin |
| Tested Applications | ELISA, WB |
| Immunogen | KLH conjugated E2 Tag peptide |
| Target | E2 tag |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Dilution Range | WB=1:1000-5000,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |