

Product Datasheet

HMG1 (phospho-Ser20/Ser24) antibody (orb4579)

Description	HMG1 (phospho-Ser20/Ser24) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HMG14 around the phosphorylation site of Ser20+Ser24 (RR(p-S)ARL(p-S)AK)
Target	phospho-HMG14 (Ser20+Ser24)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11 kDa
Uniprot ID	P05114
Expiration Date	12 months from date of receipt.