

Product Datasheet

G3BP1 (phospho-Ser232) antibody (orb4580)

Description	G3BP1 (phospho-Ser232) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human G3BP1 around the phosphorylation site of Ser232 (SS(p-S)PA)
Target	phospho-G3BP1 (Ser232)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52-60 kDa
Uniprot ID	Q13283
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.