

Product Datasheet

LMP2 antibody (Biotin) (orb458132)

Description	Rabbit polyclonal antibody against LMP2 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Proteasome 20S LMP2 (101-219/219aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23kDa
Uniprot ID	P28065
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.