

Product Datasheet

TIF1 gamma antibody (Biotin) (orb458205)

Description	TIF1 gamma antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human TIF1 gamma
Target	TIF1 gamma
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UPN9
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.