

Product Datasheet

NAT10 antibody (Biotin) (orb459224)

Description	Rabbit polyclonal antibody against NAT10 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human NAT10 (51-150/1025aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	116kDa
Uniprot ID	Q9H0A0
Dilution Range	ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.