

Product Datasheet

PTPRA (phospho-Tyr798) antibody (orb4626)

Description	PTPRA (phospho-Tyr798) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PTPRA around the phosphorylation site of Tyr798 (SD(p-Y)AN)
Target	phospho-PTPRA (Tyr798)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89 kDa
Uniprot ID	P18433
Expiration Date	12 months from date of receipt.