

Product Datasheet

Cytokeratin 19 antibody (FITC) (orb464124)

Description	Cytokeratin 19 antibody (FITC)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Cytokeratin 19
Target	Cytokeratin 19
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P08727
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.