

## Product Datasheet

### Tubulin, Alpha 4a antibody (FITC) (orb464313)

<b>Description</b>	Tubulin, Alpha 4a antibody (FITC)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC, ICC, IF
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human Tubulin, Alpha 4a
<b>Target</b>	Tubulin, Alpha 4a
<b>Preservatives</b>	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P68366</b>
<b>Dilution Range</b>	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.