

Product Datasheet

PIK3 gamma antibody (FITC) (orb465342)

Description	PIK3 gamma antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human PIK3 gamma
Target	PIK3 gamma
Preservatives	0.01M TBS(pH=7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P48736
Dilution Range	Flow-Cyt=1ug/test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.