

Product Datasheet

beta Lactoglobulin antibody (FITC) (orb466071)

Description	Rabbit polyclonal antibody against beta Lactoglobulin conjugated to FITC
Species/Host	Rabbit
Reactivity	Goat
Conjugation	FITC
Tested Applications	IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from cow Beta-lactoglobulin (11-110/178aa)
Target	LGB
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	P02754
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.