



Product Datasheet

ARHGDIB antibody (orb46680)

| | |
|----------------------------|---|
| Description | ARHGDIB antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, WB |
| Immunogen | Recombinant Human Rho GDP-dissociation inhibitor 2 protein (2-86AA) |
| Target | ARHGDIB |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P52566 |
| Dilution Range | WB:1:2000-10000, IF:1:50-200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Jurkat(lane 1),...



Immunofluorescence analysis of HepG2 cel...