

Product Datasheet

NPTX1 antibody (HRP) (orb467072)

Description	Rabbit polyclonal antibody against NPTX1 conjugated to HRP
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human NPTX1 (201-300/432aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q15818
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200
Expiration Date	12 months from date of receipt.