

Product Datasheet

Tubulin, Alpha 4a antibody (HRP) (orb474562)

| | |
|----------------------------|---|
| Description | Tubulin, Alpha 4a antibody (HRP) |
| Species/Host | Rabbit |
| Reactivity | Mouse, Rat |
| Conjugation | HRP |
| Tested Applications | ICC, IHC-Fr, IHC-P, WB |
| Immunogen | KLH conjugated synthetic peptide derived from human Tubulin, Alpha 4a |
| Target | Tubulin, Alpha 4a |
| Preservatives | Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P68366 |
| Dilution Range | WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500 |
| Expiration Date | 12 months from date of receipt. |