

Product Datasheet

HA-Tag antibody (HRP) (orb474887)

Description	HA-Tag antibody (HRP)
Species/Host	Mouse
Reactivity	All
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH conjugated HA Tag
Target	HA Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB=1:50000-1:100000
Expiration Date	12 months from date of receipt.