

Product Datasheet

EYFP-Tag antibody (HRP) (orb475129)

Description	EYFP-Tag antibody (HRP)
Species/Host	Mouse
Reactivity	All
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Recombinant EYFP protein
Target	EYFP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P21578
Dilution Range	WB=1:1000-5000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.