



## Product Datasheet

### TNFAIP3 antibody (orb47519)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | TNFAIP3 antibody   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, IF, WB  |
| <b>Immunogen</b>           | Recombinant Human Tumor necrosis factor alpha-induced protein 3 protein (97-183AA) |
| <b>Target</b>              | TNFAIP3  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                        |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                      |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <b>P21580</b>  |
| <b>Dilution Range</b>      | WB:1:2000-10000, IF:1:50-500   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Western blot analysis of RAW264.7(lane 1...



Immunofluorescence analysis of HeLa cell...