



Product Datasheet

ATG4B antibody (orb47863)


| | |
|----------------------------|---|
| Description | ATG4B antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IHC, WB |
| Immunogen | Recombinant Human Cysteine protease ATG4B protein (350-400AA) |
| Target | ATG4B |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9Y4P1 |
| Dilution Range | WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of human co...



Western blot analysis of Mouse brain tis...



Immunohistochemical staining of human pa...