

Product Datasheet

HP1 gamma (phospho-Ser83) antibody (orb4790)

Description	HP1 gamma (phospho-Ser83) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HP1 gamma around the phosphorylation site of Ser83 (RK(p-S)LS)
Target	Phospho-HP1 gamma (Ser93)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	Q13185
Expiration Date	12 months from date of receipt.