

Product Datasheet

Smad1 (phospho-Ser463/465) antibody (HRP) (orb480042)

| | |
|----------------------------|--|
| Description | Smad1 (phospho-Ser463/465) antibody (HRP) |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | HRP |
| Tested Applications | ELISA, ICC, IHC-Fr, IHC-P, WB |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human Smad1/5 around the phosphorylation site of Ser463/465 |
| Target | Phospho-Smad1/5 (Ser463 + Ser465) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q15797 |
| Dilution Range | WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |