

Product Datasheet

PRE7 antibody (HRP) (orb480287)

Description	Rabbit polyclonal antibody against PRE7 conjugated to HRP
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PRS3 (41-140/318aa)
Target	PRPS1L1
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35 kDa
Uniprot ID	P21108
Dilution Range	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.