


Product Datasheet

C-Myc Tag antibody (orb48094)

Description	Mouse monoclonal to CMyc Tag.c Myc tag antibody is part of the Tag series of antibodies, the best qu...
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:3000-1:10000, IF:1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK 293 cell ly...