



Product Datasheet

EDNRB antibody (orb48152)


Description	EDNRB antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant Human Endothelin receptor type B protein (27-101AA)
Target	Tissue Specificity: Expressed in placental stem villi vessels, but not in cultured placental villi smooth muscle cells.
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 95%, Protein G purified
MW	Observed band size: 50 kDa
Uniprot ID	P24530
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart(lan...



Immunofluorescence analysis of PC-3 cell...



Immunohistochemical staining of human lu...