



Product Datasheet

GARS antibody (orb48232)


| | |
|----------------------------|---|
| Description | GARS antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IHC, IP, WB |
| Immunogen | Recombinant Human Glycine--tRNA ligase protein (1-389AA) |
| Target | GARS |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P41250 |
| Dilution Range | WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200, IP:1:200-2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of THP-1(lane 1), ...



Immunohistochemical staining of human te...



Immunofluorescence analysis of HepG2 cel...