

Product Datasheet

BMP11 antibody (FITC) (orb48242)

Description	BMP11 antibody (FITC)
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	FITC
Immunogen	Recombinant Pig Growth differentiation factor 11 protein (296-404AA)
Target	GDF11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BG54
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.