

Product Datasheet

Bcr (phospho-Tyr177) antibody (orb4827)

Description	Bcr (phospho-Tyr177) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesisedphosphopeptide derived from human Bcr around the phosphorylation site of Tyr177 (PF(p-Y)VN)
Target	Phospho-Bcr (Tyr177)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	143 kDa
Uniprot ID	P11274
Expiration Date	12 months from date of receipt.