

Product Datasheet

Cellulase antibody (HRP) (orb482718)

Description	Rabbit polyclonal antibody to Cellulase
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	Cellulase from Trichoderma reesei ATCC 26921
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35kDa
CAS Number	9012-54-8
Dilution Range	ELISA: 1:100-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.