



Product Datasheet

BTK antibody (orb48294)


Description	Rabbit polyclonal antibody to BTK
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human BTK (NP_000052.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q06187
Entrez	695
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



IHC-P staining of rat lung tissue using ...



Immunohistochemical analysis of formalin...