



Product Datasheet

Galectin 3 antibody (orb48437)


Description	Rabbit polyclonal antibody to Galectin 3
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Galectin 3/Galectin 3/LGALS3(NP_002297.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	P17931
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-50
Expiration Date	12 months from date of receipt.



IHC-P staining of mouse kidney tissue us...



Immunofluorescence analysis of MCF7 cell...



Immunohistochemistry of paraffin-embedde...