

Product Datasheet

SULT1A1 antibody (orb48551)

Description Rabbit polyclonal to SULT1A1

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant fusion protein of human SULT1A1(NP_803565.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 34kDa


Uniprot ID [P50225](#)

Dilution Range WB: 1:500-2000, IHC-P: 1:50-200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of mouse kidney



Immunohistochemical staining of rat brain



Immunohistochemical staining of rat kidney