



Product Datasheet

SIGLEC9 antibody (orb4859)


Description	SIGLEC9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human CD329 (201-300/499 aa)
Target	SIGLEC8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q9NYZ4
Dilution Range	Flow-Cyt=1µg /test, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Raji cell us...



Flow cytometric analysis of Raji cell us...



Immunohistochemical staining of human lu...