

Product Datasheet

EGFR antibody (PE) (orb485910)

Description	EGFR antibody (PE)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from human EGFR
Target	EGFR
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00533
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.