

## Product Datasheet

### GBAS antibody (PE) (orb486790)

<b>Description</b>	Rabbit polyclonal antibody against GBAS conjugated to PE
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	PE
<b>Tested Applications</b>	ICC, IF
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human GBAS
<b>Target</b>	NIPSNAP2
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300.
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	- 100 µg (in 100 µl): 1 mg/ml
<b>Storage</b>	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	34kDa
<b>Uniprot ID</b>	<a href="#">O75323</a>
<b>Dilution Range</b>	IF=1:100-500,ICC=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.