



## Product Datasheet

### IL17A antibody (orb48920)

<b>Description</b>	Rabbit polyclonal to IL17A.IL-17A is a cystine-linked homodimeric pro-inflammatory cytokine produced...
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant fusion protein of human IL17A(NP_002181.1).
<b>Preservatives</b>	PBS with 0.05% proclin300,50% glycerol,pH7.3.
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	18-25KDa
<b>Uniprot ID</b>	<a href="#">Q16552</a>
<b>Dilution Range</b>	WB 1:500 - 1:2000, IF/ 1:50 - 1:200
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Western blot analysis of THP-1 cells usi...