



Product Datasheet

TARDBP antibody (orb49083)


Description	Rabbit polyclonal antibody to TARDBP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human TDP-43/TARDB(NP_031401.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42kDa
Uniprot ID	Q13148
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:20-100, IP: 1:20-50R, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of HeLa cell...



Immunohistochemistry of paraffin-embedde...