

Product Datasheet

Human ACE2 protein (orb49193)

Description	Recombinant Human Angiotensin-Converting Enzyme 2 is produced by our Mammalian expression system and...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mMTris,300mMNacl,1mMZncl2,10%Glycerol,PH7.4
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	Q9BYF1
Expiration Date	6 months from date of receipt.