

Product Datasheet

Human APE protein (orb49206)

Description	Recombinant Human Apurinic-Apyrimidinic Endonuclease 1 is produced by our E.coli expression system a...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 10mM HEPES, 100mM KCl, 50% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Uniprot ID	AAH02338.1
Expiration Date	6 months from date of receipt.