

Product Datasheet

Human BMP2 protein (orb49220)

Description	Recombinant Human/Mouse/Rat Bone Morphogenetic Protein 2 is produced by our E.coli expression system...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 10mM HAc-NH ₄ Ac, 4% D-Mannitol, pH 4.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	P12643
Expiration Date	6 months from date of receipt.