

Product Datasheet

JunD (phospho-Ser255) antibody (orb4925)

Description	JunD (phospho-Ser255) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human JunD around the phosphorylation site of Ser255 (GE(p-S)PP)
Target	phospho-JunD (Ser255)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	P17535
Expiration Date	12 months from date of receipt.