

Product Datasheet

Scavenger Receptor BI + BII antibody (PE) (orb492991)

Description	Rabbit polyclonal antibody against Scavenger Receptor BI + BII conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Scavenger Receptor BI + BII (31-130/552aa)
Target	SCARB1
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61kDa
Uniprot ID	Q8WTV0
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.