

Product Datasheet

Human Cyclophilin C protein (orb49331)

Description	Recombinant Human Peptidyl-Prolyl cis-trans Isomerase C is produced by our E.coli expression system ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-Trx&His
Note	For research use only
Application notes	Trx, 6His tag at the N-terminus.
Uniprot ID	P45877
Expiration Date	6 months from date of receipt.