

Product Datasheet

Human GOT1 protein (orb49408)

Description	Recombinant Human Glutamate Oxaloacetate Transaminase 1 is produced by our E.coli expression system ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 2mM DTT, 20% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	P17174
Expiration Date	6 months from date of receipt.