

Product Datasheet

Human IL1F10 protein (orb49459)

Description	Recombinant Human Interleukin-1 Family Member 10 is produced by our E.coli expression system and the...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 6.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	AAI03967.1
Expiration Date	6 months from date of receipt.