

Product Datasheet

Human IL11 protein (orb49463)

Description	Recombinant Human Interleukin-11 is produced by our Yeast expression system and the target gene enco...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 2% Glycine, pH 7.2.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	P20809
Expiration Date	6 months from date of receipt.