

Product Datasheet

Human Kallikrein 11 protein (orb49488)

Description	Recombinant Human Kallikrein 11 is produced by our Mammalian expression system and the target gene e...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 2mM CaCl ₂ , pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	AAH22068
Expiration Date	6 months from date of receipt.