

Product Datasheet

BRMS1 antibody (PE) (orb495094)

Description	Rabbit polyclonal antibody against BRMS1 conjugated to PE
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human BRMS-1 (151-246/246aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	Q9HCU9
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.