

Product Datasheet

Human PBLD protein (orb49581)

Description	Recombinant Human Phenazine Biosynthesis-Like Domain-Containing Protein is produced by our E.coli ex...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 30% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Uniprot ID	P30039
Expiration Date	6 months from date of receipt.