
Product Datasheet

RNA polymerase II CTD repeat YSPTSPS (phospho-Ser5) antibody (PE) (orb496617)

Evolution Bioproducts.

RNA polymerase II CTD repeat

Description	RNA polymerase II CTD repeat (PE)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RNA polymerase II CTD around the phosphorylation site of Ser5
Target	RNA polymerase II CTD repeat YSPTSPS (phospho S5)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24928
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.