

Product Datasheet

Human HR6A protein (orb49713)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 A is produced by our E.coli expression system and ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-GST&His
Note	For research use only
Application notes	GST, 6His tag at the N-terminus.
Uniprot ID	P49459
Expiration Date	6 months from date of receipt.